**steriset**<sup>®</sup>


**Sterilcontainer**

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CE

**Consequently safe**

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WAGNER 

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Honored\* as "Exemplary Product for Safety in Medicine" in the "Golden SMM '98" competition

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\* SMM = German Society for "Safe Materials in Medicine e.V."

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### Proven World-Wide

WAGNER is a **leading supplier of Sterilization Containers** for hospitals.

The **SteriSet System** offers the results of 50 years of marketing experience and consistent fulfillment of the requirements in clinical practice for providing modern sterilized goods: A **high-quality, multiple patented, sterile container system** that solves the essential requirements such as safety, simple handling, economy

and environmental compatibility in a multi-patented manner.

### Quality

SteriSet is manufactured in Germany by trained and qualified personnel under a TÜV certified quality management system according to DIN EN ISO 9001/46001.

### Safety

SteriSet was designed according to the principle of passive safety: Less or improved handling to reduce the "human" error factor as well as more protection from incalculable environmental influences.

### Flexibility

SteriSet is offered in three sizes, half, three-quarter and full size containers. In addition, there are five different heights available, sealable or self-locking ThermoLoc closure, a condensate drain, and color-coded cover and handle to differentiate each container.

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## Perfect Protection

Even the best bio-barriers are of no value when they are subjected to influences that can be detrimental to the barrier characteristics – and thereby the maintenance of sterility – during sterilization, storage or transport.

**SteriSet comes standard with a protective cover** that protects the actual germ barrier effectively against detrimental environmental influences.

Although such a container is nearly “impermeable” in the closed state, it can be steam-sterilized with the cover in place and sealed.

This “Protective Cover” design improves storage and transport capabilities and thereby provides extraordinarily high hygienic safety.

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


During steam sterilization (pressure difference), steam-air exchange can take place without problems through the circumferential gap between the cover and the container.

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With the container open, the inner barrier cover can be removed from the outer cover at any time (e.g. for checking or cleaning).

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## **Aseptic Presentation by Design**

Even with proper storage and correct transporting, the packaging of sterile goods becomes more and more contaminated on the outside, especially the horizontal surfaces – depending on the quality and duration of storage.

The increased handling risk is already drastically reduced by the protective design of the SteriSet containers, but how can they be safely opened?

**The optional “SteriSet S-Model”** provides a 100% consistent solution for **organized** aseptic presentation since **opening takes place in two steps:**

■ After unsealing and removal of the outer protective cover, the inner lid (= the actual sterile barrier) still remains closed by the internally located additional closure.

Since the barrier cover was protected against sedimentation during the entire storage time, it is “not contaminated”, i.e. with only approx. 1% of external contamination in comparison with the outer cover or, for example, with a filter that is not protected against sedimentation.

■ It can then be opened as the **second stage**, e.g., within the operating room and without danger of contaminating the environment or the contents: **Organized aseptic presentation by design!**

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## Really Locked!

What is the use of a “seal” on a closure when the complete container including the barrier and the gasket is not protected against undetected opening?

SteriSet containers, with the patented **ThermoLoc closure**, provide a solution to this required sealing in an outstandingly elegant manner: **They “seal” themselves during sterilization**. By means of the temperature effect, a color-marked catch is slid in front of the closure latch such that swinging it open is not possible. This “sealed” state is visually recognizable – and provides above all a **perceptible blockage of the closure**: The barrier cover together with the gasket is then securely “locked” with the lower part – and reusable without handling and without consumables.

The container is opened in **two** steps:

- 1 Intentional unsealing of the closure blockage,
- 2 Swinging open of the closure latches.

This eliminates the meticulous visual checking of both container fronts required by other systems for intactness of the sealing, which is an invaluable advantage in the often hectic operating room routine.

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
The latch of the outer cover is securely blocked against being swung open by the stable, colored catch of the ThermoLoc closure.

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The ThermoLoc closure in sealed and unsealed condition.

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ThermoLoc containers are unsealed after cooling by a simple hand motion (upward pressing) on the colored catch. In this manner, the closure is also automatically prepared for the next sealing (sterilization).

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## Wear-free Bio-barrier...

SteriSet **Valve** Containers are containers with reusable, **wear-free** pressure / vacuum valves of electrolytically polished **stainless steel**.

During steam sterilization, the valves automatically follow the pressure changes that occur and enable the unrestricted flow of the steam / air in and out of the container.

After removal from the sterilizer, the valves are closed by spring force and offer an effective micro-biological barrier (tortuous path concept).

The SteriSet **Valve** System is proven technology that has been used for several decades. It **eliminates annoying filter changing or checking**, reduces the handling effort, is really long-lived (stainless steel: 10 year guarantee!) and thereby **extremely cost efficient**.

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



**The SteriSet valves made of polished stainless steel are long-lived, wear-free and self-cleaning (steam-jet principle). With the outer lid removed, they can be quickly checked for their functionality.**

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## ... Or Filter After All?

As an alternative to the valve version is the SteriSet **Filter** System.

With conventional (flexible, thin) filter holders there is always a risk that, – with a slight deformation only of the holder lid (see illustration) occurring before, during or after sterilization, – the filter member might not correctly seal to the lid – which would result in the loss of the sterile barrier. The **filter version of the SteriSet System** therefore deliberately places great value on:

- Highest possible safety by means of a centrally attached, spring-suspended, **extremely robust filter holder**,
- secure, sealed pressing of the filter against the inner cover (unevenness of up to 5 mm can be compensated),
- quick replacing of the filter with one-hand operation.

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**Danger** with already only lightly deformed conventional filter holders (a) or covers (b): with a filter not correctly pressed on, there is no sterile barrier!

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**Advantage:** only one size (round) filter, suitable for all containers, fast and easy to exchange and control.

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## Safe Drying!

Experience indicates that, with increasing weight, the probability increases that residual moisture will remain in the container after sterilization.

**The SteriSet containers with condensate drain remove excess condensate** – especially under unfavorable conditions.

For example, if a load of 10 kg is exceeded (orthopedics, heavy metal or plastic load) or if unfavorable sterilization conditions exist (wet steam, older sterilization equipment...), a container with closed bottom quickly reaches its limits during drying.

With the **patented condensate drain**, the probability of drying even heavy metal loads is clearly improved: This is achieved by a thermally controlled bottom drain – placed at the deepest point of the base tray – that opens **during**

**sterilization** when 130 °C are reached (as an alternative, a 121 °C version is also available) and thereby **excess condensate in fluid form is removed from the container** during the sterilization process.

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
**The patented condensate valve in closed and opened state (functional diagram).**

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

Still during the vacuum drying phase, the valve of the bottom drain closes **the first time the temperature goes below 110 °C** with high contact pressure and seals tightly.

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The encapsulated, one-piece condensate valve can be removed and replaced with a simple press/turn motion.

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## Expandable

A system is complete only when suitable accessories are also available.

Supplementary to the containers, the SteriSet system provides you with the **SteriSet Small-Set program** as well as a comprehensive **accessory program** that ranges from consumables, **strainer pans** and divider elements to **tables and cars**, everything from one source.

We would be pleased to send you a detailed product catalog for the individual product lines.

In addition, all catalogs are available to you at [www.wagner-steriset.de](http://www.wagner-steriset.de) in PDF versions for downloading.

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 **Leader in the Supplying of Sterile Goods for more than 50 Years**

WAGNER is a medium-sized enterprise in the medical equipment sector with headquarters in Munich, Germany. Originally a craftsman's workshop, today the worldwide active WAGNER GmbH with its traditionally innovative products has developed to become a **leading supplier of sterile containers** for hospitals.

The consistent application of our strategy "Quality in manufacturing, innovation in design as well as cooperation in the markets" enables us today, **in teamwork with a worldwide network of distributors**, to offer a high quality, multiple-patented container system in various versions.

**SteriSet containers fulfil the requirements of the relevant DIN/EN standards** of the "Medicine Product Laws" (CE certification) as well as the entry requirements for the US market (FDA 510K Premarketing Notification).

The products are manufactured in Germany by trained and qualified personnel under a TUV certified quality management system acc. to DIN EN ISO 9001/46001 and are being constantly refined.

**For further current information** as well as all product catalogs, operating instructions, answers to frequently asked questions, etc. visit our site on the Internet under [www.wagner-steriset.de](http://www.wagner-steriset.de)

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SteriSet® Containers are manufactured in Germany, covered by a TÜV-certified Quality-Management-System acc. DIN EN ISO 9001 and in addition DIN EN 46001 (for medical devices).

The design is protected by multiple international patents.

**SteriSet® Containers meet the following requirements of the:**

- ◆ European Medical Device Directive (CE-mark declared)
- ◆ European Standards EN 868-1/EN 868-8
- ◆ German Standards (Former DIN 58952/1 and actual DIN 58953/9)
- ◆ FDA 510K Premarketing Notification for USA

SteriSet® and ThermoLoc® are registered trademarks

Technical modifications reserved

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© 2003
**Wagner GmbH**